



## 30 MINUTE FIRE RATED HAZARDOUS STORAGE CABINET

### SL-B802640A



Store LAB (Type 30) with Fire Resistance >30 minutes Flammable/Corrosive/Acid/Toxic/Pesticide Cabinet, 4 Gallon, 15 Litre, Self Closing, 1 Door, 1 Shelf, Grey Body RAL7035, Daffodil Yellow Door RAL1007, Levering Feet, Auto Thermal Closing Mechanism, Height-adjustable shelves at intervals of 32 mm, Levering Feet for levelling, With Base Sump Tray and Perforated Cover Plate at Base, Cylinder Lock on Door, Door Hinge to Right. Markings according to DIN EN 14470-1 / TRGS 510 / ISO 3864 / ASR A1.3.

**Specifications:** 4G / 15L - Packed Dimensions (HxWxD) 800 x 1200 x 800mm



SPECIFICATIONS		
Item Code	External Dimensions	Weight kg
SL-B802640A	600mm(H) x 600mm(W) x 570mm(D)	63 kg
	Internal Dimensions	Packing Weight kg
	425mm(H) x 472mm(W)x 470mm(D)	80 kg

### SL-B802609A



Store LAB (Type 30) with Fire Resistance >30 minutes Flammable/Corrosive/Acid/Toxic/Pesticide Cabinet, 12 Gallon, 45 Litre, Self Closing, 2 Doors, 1 Shelf, Grey Body RAL7035, Daffodil Yellow Door RAL1007, Levering Feet, Auto Thermal Closing Mechanism, Height-adjustable shelves at intervals of 32 mm, Levering Feet for levelling, With Base Sump Tray and Perforated Cover Plate at Base, Cylinder Lock on Doors, Door Hinge to Right. Markings according to DIN EN 14470-1 / TRGS 510 / ISO 3864 / ASR A1.3

**Specifications:** 12G / 45L - Packed Dimensions (HxWxD) 800 x 1200 x 800 mm



SPECIFICATIONS		
Item Code	External Dimensions	Weight kg
SL-B802609A	600mm(H) x 1100mm(W) x 570mm(D)	110 kg
	Internal Dimensions	Packing Weight kg
	425mm(H) x 986mm(W)x 455mm(D)	125 kg

**SL-B802608A**

Store LAB (Type 30) with Fire Resistance >30 minutes Flammable/Corrosive/Acid/ Toxic/Pesticide Cabinet, 15 Gallon, 57 Litre, Self Closing, 1 Door, 2 Shelves, Grey Body RAL7035, Daffodil Yellow Door RAL1007, Levering Feet, Auto Thermal Closing Mechanism, Height-adjustable shelves at intervals of 32 mm, Levering Feet for levelling, With Base Sump Tray and Perforated Cover Plate at Base, Cylinder Lock on Door, Door Hinge to Right. Markings according to DIN EN 14470-1 / TRGS 510 / ISO 3864 / ASR A1.3.

**Specifications:** 15G / 57L - Packed Dimensions (HxWxD) 1500 x 1200 x 800 mm

**SPECIFICATIONS**

Item Code	External Dimensions	Weight kg
<b>SL-B802608A</b>	1315mm(H) x 595mm(W) x 595mm(D)	122 kg
	Internal Dimensions	Packing Weight kg
	1108mm(H) x 475mm(W)x 515mm(D)	140 kg

**SL-B802603A**

Store LAB (Type 30) with Fire Resistance >30 minutes Flammable/Corrosive/Acid/ Toxic/Pesticide Cabinet, 27 Gallon, 103 Litre, Self Closing, 1 Door, 3 Shelves, Grey Body RAL7035, Daffodil Yellow Door RAL1007, Levering Feet, Auto Thermal Closing Mechanism, Height-adjustable shelves at intervals of 32 mm, Levering Feet for levelling, With Base Sump Tray and Perforated Cover Plate at Base, Cylinder Lock on Door, Door Hinge to Right. Markings according to DIN EN 14470-1 / TRGS 510 / ISO 3864 / ASR A1.3.

**Specifications:** 27G / 103L - Packed Dimensions (HxWxD) 2200 x 1200 x 800 mm

**SPECIFICATIONS**

Item Code	External Dimensions	Weight kg
<b>SL-B802603A</b>	1935mm(H) x 595mm(W) x 595mm(D)	147 kg
	Internal Dimensions	Packing Weight kg
	1600mm(H) x 475mm(W)x 515mm(D)	165 kg

**SL-B802604A**

Store LAB (Type 30) with Fire Resistance >30 minutes Flammable/Corrosive/Acid/ Toxic/Pesticide Cabinet, 60 Gallon, 227 Litre, Self Closing, 2 Door, 3 Shelves, Grey Body RAL7035, Daffodil Yellow Door RAL1007, Levering Feet, Auto Thermal Closing Mechanism, Height-adjustable shelves at intervals of 32 mm, Levering Feet for levelling, With Base Sump Tray and Perforated Cover Plate at Base, Cylinder Lock on Doors, Door Hinge to Right. Markings according to DIN EN 14470-1 / TRGS 510 / ISO 3864 / ASR A1.3.

**Specifications:** 60G / 227L - Packed Dimensions (HxWxD) 2200 x 1200 x 800 mm

**SPECIFICATIONS**

Item Code	External Dimensions	Weight kg
<b>SL-B802604A</b>	1935mm(H) x 1195mm(W) x 595mm(D)	225 kg
	Internal Dimensions	Packing Weight kg
	1600mm(H) x 1075mm(W)x 515mm(D)	240 kg