

## Waste Bags

Our Range of Contaminated/Hazardous Waste Bags, Clinical Waste bags & Bio-Hazard Bags; All bags are manufactured to industry requirements for thickness and disposal.

Made from LDPE with good tear and puncture resistance, they are suitable for the disposal of biologically hazardous waste. Comes with highly visible "Hazardous / Biohazard / Clinical Waste" labels and cable ties for easy sealing and disposal.

### SPECIFICATIONS

Item	Details
X-BAT	HAZARDOUS WASTE BAG 600 X 700 X 0.04MM THICKNESS - LDPE WITH CABLE TIES
X-CWB	CLINICAL WASTE BAG 700 X 550mm X 0.05MM THICKNESS - LDPE WITH CABLE TIES
X-BWB120	BIOHAZARD WASTE BAG YELLOW - TO SUIT 120L/140L TWO WHEEL BIN-(700+300)MM x1200MM - 0.05MM-LDPE
X-BWB240	BIOHAZARD WASTE BAG YELLOW - TO SUIT 240L TWO WHEEL BIN - (970 +300)MM x 1500MM - 0.05MM-LDPE

